



Quality Medical Products

Zefon Arterial Anchor Bandage



Secures Arterial Lines in One Simple, Sterile, & Suture-free Procedure

LateX
Free

Unlike conventional "tape & gauze" securement techniques, the Zefon Arterial Line Securement Device, with its Clear-vue Dome, protects the catheter entry point from direct trauma and provides for continual observation of the site, allowing for quicker detection and correction of potential problems. The E-Z Lok Channel safely secures the catheter in a proper fixed position and eliminates the need for suturing, thus reducing the risk of infection and improving patient comfort.

- Clear-vue dome accommodates cannulas and pressure lines with long hubs as well as arterial catheters with "wing sets"
- Small dome footprint provides excellent fit on patients with small arms/wrists
- Strong adhesive bandage holds but won't harm or irritate patients skin

Download Zefon Arterial Line Securement Device Brochure

Features and Benefits

- Clear-Vue Dome provides physical protections & direct site visibility
- Vents provide for air circulation over the site, enhancing the healing process.
- E-Z Lok Channel securely holds catheter at a proper angle.
- Adhesive backing securely holds bandage in place without irritating patients skin.
- Absorbent Pad draws excess fluid away from puncture site.

- Quick-clips hold extension line neatly in place.



Dr. Z Says.....

The Zefon arterial line holder is designed to hold the catheter at the proper angle and in a fixed position, eliminating the need for suturing, reducing risk of infection and improving patient comfort.

Ordering Information

Product #	Description
ZAB-4100	Zefon Arterial Line Holder, Sterile, Individually Packaged, 50/bx

All Zefon medical products are sold exclusively through specialty distributors.

Product availability may vary by country/region. Please contact your local Zefon distributor for product, pricing and availability information in your area. A list of distributors can be found in the **Where to Buy** section on this website.

Accessories & Related Products



GRIP-LOK™ Universal Catheter and Tubing Securement



GRIP-LOK™ H Hydrocolloid Universal Catheter and Tubing Securement

IMA
Independent Medical Associates

www.I-MA.com
info@I-MA.com
p: 888.548.4462
f: 888-548-1462