

ComfortCheck™



Disposable Blood Pressure Cuffs

For over ten years Zefon International has been designing and manufacturing superior blood pressure cuffs for leading blood pressure measurement instrument manufacturers. With a complete range of sizes, materials and connectors, Zefon brand cuffs represent the best value to bring consistent performance, patient comfort and economy to your medical setting.

Features

- Disposable; Eliminates costly disinfection procedures and contamination concerns.
- Cost Effective; Economical enough for single patient use.
- Increase infection control with single patient use cuffs.
- Soft white fabric material; comfortable for patients.
- Meets AAMI and AHA Cuff Sizing Recommendations.
- Color coded printing on each cuff makes it easy to identify cuff size.
- Latex Free.
- Designed as a direct replacement for OEM cuffs provided by monitor manufacturers.
- Full range of connectors and adapters allow standardization of cuffs across the whole hospital.
- Easy connection to electronic NIBP monitors or manual blood pressure measurement systems.

Latex Free **CE**




Zefon International Inc.
5350 SW 1st Lane, Ocala, FL 34474 USA
Tel: +1 3528548080; Fax: +1 3528547480
www.zefon.com

MM00001 Rev 5 2008-06

EC REP Blue Box Medical Limited
Unit 29, New Forest Enterprise Centre,
Rushington Business Park; Chapel Lane
Totton, Southampton, SO40 9LA United Kingdom
Phone: +023 8066 9000; Fax: +023 8066 9002

Ordering Information

Single Tube Cuffs




		Monitor				
		CAS 740/750 Criticare 508 Datascopes Accutor Plus & Passport 2 Philips/HP Siemens Spacelabs	Colin Datascopes Passport XG & Expert Nihon Kodon Lifescopes 9 & 12 Welch Allyn Vital Signs/Atlas (Before July 2000)	Welch Allyn Vital Signs/Atlas (after July 2000)	Welch Allyn Protocol Propaq	
Cuff Size	Color Code	Limb Circumference	1-Tube Bayonet Connector	1-Tube Female Luer Connector	1-Tube Screw Type Connector	1-Tube Female Twist Lock Connector
Child	Green	11-22cm	ZCC4401HP	ZCC4401F	ZCC4401D	ZCC4401MQF
Small Adult	Light Blue	16-28cm	ZCC4301HP	ZCC4301F	ZCC4301D	ZCC4301MQF
Adult	Navy Blue	26-33cm	ZCC4201HP	ZCC4201F	ZCC4201D	ZCC4201MQF
Adult Long*	Navy Blue	26-33cm	ZCC4221HP	ZCC4221F	ZCC4221D	ZCC4221MQF
Large Adult	Red	33-41cm	ZCC4101HP	ZCC4101F	ZCC4101D	ZCC4101MQF
Large Adult Long*	Red	33-41cm	ZCC4121HP	ZCC4121F	ZCC4121D	ZCC4121MQF
Thigh	Brown	39-55cm	ZCC4001HP	ZCC4001F	ZCC4001D	ZCC4001MQF

* Available late summer 2008

All cuffs are packaged with 10 cuffs per box

Double Tube Cuffs



		Monitor				Manual
		Critikon Dinamaf Criticare 507 Invivo Research Nihon Kodon (Lifescopes L) Spacelabs	Datex-Ohmeda	GE Medical Eagle 4000 Tram 100/200/300		
Cuff Size	Color Code	Limb Circumference	2-Tube Screw Type Connector	2-Tube Male/Female Twist Lock Connector	2-Tube Female Twist Lock Connector	2-Tube Manual Inflation with Male Luer
Child	Green	11-22cm	ZCC4402D	ZCC4402MQMF	ZCC4402MQF	ZCC4412L
Small Adult	Light Blue	16-28cm	ZCC4302D	ZCC4302MQMF	ZCC4302MQF	ZCC4312L
Adult	Navy Blue	26-33cm	ZCC4202D	ZCC4202MQMF	ZCC4202MQF	ZCC4212L
Adult Long*	Navy Blue	26-33cm	ZCC4222D	ZCC4222MQMF	ZCC4222MQF	ZCC4232L
Large Adult	Red	33-41cm	ZCC4102D	ZCC4102MQMF	ZCC4102MQF	ZCC4112L
Large Adult Long*	Red	33-41cm	ZCC4122D	ZCC4122MQMF	ZCC4122MQF	ZCC4132L
Thigh	Brown	39-55cm	ZCC4002D	ZCC4002MQMF	ZCC4002MQF	ZCC4012L

* Available late summer 2008

All cuffs are packaged with 10 cuffs per box

Distributed By:



5350 SW 1st Lane,
Ocala, FL 34474 USA
Tel: +1 3528548080
www.zefon.com

IMA

Independent Medical Associates

11733 66th St N • #113 • Largo, FL • 33773 • P: 888-548-4462

• F: 888-548-1462 • info@I-MA.com • www.I-MA.com