

DuraCheck™

Reusable Bladderless Blood Pressure Cuffs



For over ten years Zefon International has been designing and manufacturing superior blood pressure cuffs for leading blood pressure measurement instrument manufacturers. With a complete range of sizes, materials and connectors, Zefon brand cuffs represent the best value to bring consistent performance, patient comfort and economy to your medical setting.

Features

- Reusable; Durable construction for long life and multiple patient use.
- Bladderless design is accurate and easy to clean.
- Made from durable nylon material that is comfortable and strong.
- Meets AAMI and AHA Cuff Sizing Recommendations.
- Color coding makes it easy to identify cuff size.
- Latex Free.
- Complete selection of connectors and adapters allow standardization of cuffs across the whole hospital.
- Easy connection to electronic NIBP monitors or manual blood pressure measurement systems.

Latex Free CE




5350 SW 1st Lane, Ocala, FL 34474 USA
Tel: +1 3528548080; Fax: +1 3528547480
www.zefon.com

MM00002 Rev 5 2008-06

EC REP Blue Box Medical Limited
Unit 29, New Forest Enterprise Centre,
Rushington Business Park; Chapel Lane
Totton, Southampton, SO40 9LA United Kingdom
Phone: +023 8066 9000; Fax: +023 8066 9002

Ordering Information

Single Tube Cuffs




			Monitor			
			CAS 740/750 Criticare 508 Datascope Accutor Plus & Passport 2 Philips/HP Siemens Spacelabs	Colin Datascope Passport XG & Expert Nihon Kodan Lifescope 9 & 12 Welch Allyn Vital Signs/Atlas (Before July 2000)	Welch Allyn Vital Signs/Atlas (after July 2000)	Welch Allyn Protocol Propaq
Cuff Size	Color Code	Limb Circumference	1-Tube Bayonet Connector	1-Tube Female Luer Connector	1-Tube Screw Type Connector	1-Tube Female Twist Lock Connector
Child	Green	11-22cm	ZDC4401HP	ZDC4401F	ZDC4401D	ZDC4401MQF
Small Adult	Light Blue	16-28cm	ZDC4301HP	ZDC4301F	ZDC4301D	ZDC4301MQF
Adult	Navy Blue	26-33cm	ZDC4201HP	ZDC4201F	ZDC4201D	ZDC4201MQF
Adult Long*	Navy Blue	26-33cm	ZDC4221HP	ZDC4221F	ZDC4221D	ZDC4221MQF
Large Adult	Red	33-41cm	ZDC4101HP	ZDC4101F	ZDC4101D	ZDC4101MQF
Large Adult Long*	Red	33-41cm	ZDC4121HP	ZDC4121F	ZDC4121D	ZDC4121MQF
Thigh	Dark Blue	39-55cm	ZDC4001HP	ZDC4001F	ZDC4001D	ZDC4001MQF

* Available late summer 2008

All cuffs are packaged with 5 cuffs per box

Double Tube Cuffs



			Monitor			
			Critikon Dinamap Criticare 507 Invivo Research Nihon Kodan (Lifescope L) Spacelabs	Datex-Ohmeda	GE Medical Eagle 4000 Tram 100/200/300	Manual
Cuff Size	Color Code	Limb Circumference	2-Tube Screw Type Connector	2-Tube Male/Female Twist Lock Connector	2-Tube Female Twist Lock Connector	2-Tube Manual Inflation with Male Luer
Child	Green	11-22cm	ZDC4402D	ZDC4402MQMF	ZDC4402MQF	ZDC4412L
Small Adult	Light Blue	16-28cm	ZDC4302D	ZDC4302MQMF	ZDC4302MQF	ZDC4312L
Adult	Navy Blue	26-33cm	ZDC4202D	ZDC4202MQMF	ZDC4202MQF	ZDC4212L
Adult Long*	Navy Blue	26-33cm	ZDC4222D	ZDC4222MQMF	ZDC4222MQF	ZDC4232L
Large Adult	Red	33-41cm	ZDC4102D	ZDC4102MQMF	ZDC4102MQF	ZDC4112L
Large Adult Long*	Red	33-41cm	ZDC4122D	ZDC4122MQMF	ZDC4122MQF	ZDC4132L
Thigh	Dark Blue	39-55cm	ZDC4002D	ZDC4002MQMF	ZDC4002MQF	ZDC4012L

* Available late summer 2008

All cuffs are packaged with 5 cuffs per box



5350 SW 1st Lane,
Ocala, FL 34474 USA
Tel: +1 3528548080
www.zefon.com

Distributed By:



www.I-MA.com
info@I-MA.com
p: 888.548.4462
f: 888-548-1462