

# Zefon Gas Sampling Lines

Economical  
Disposable, Single Patient Use  
Non-Sterile  
Sizes Available for Any Monitor  
100% Flow Tested  
Top Quality Medical Grade Tubing  
Filters & Nafion® Dryer Tube Options  
Custom Configurations Available



**LATEX FREE**

## Ordering Information

### 8' PVC Gas Sampling Lines

Product #	Tube ID	Connectors
ZML-9300	.040"	Male/Male
ZML-9400	.040"	Male/Female
ZML-9330	.050"	Male/Male
ZML-9430	.050"	Male/Female
ZML-9320	.060"	Male/Male
ZML-9420	.060"	Male/Female

### 10' PVC Gas Sampling Lines

Product #	Tube ID	Connectors
ZML-9500	.040"	Male/Male
ZML-9600	.040"	Male/Female
ZML-9530	.050"	Male/Male
ZML-9630	.050"	Male/Female
ZML-9520	.060"	Male/Male
ZML-9620	.060"	Male/Female

### 8' Co-Extruded PE/PVC Gas Sampling Lines

Product #	Tube ID	Connectors
ZML-9350	.040"	Male/Male
ZML-9450	.040"	Male/Female
ZML-9340	.047"	Male/Male
ZML-9440	.047"	Male/Female

### 10' Co-Extruded PE/PVC Gas Sampling Lines

Product #	Tube ID	Connectors
ZML-9550	.040"	Male/Male
ZML-9650	.040"	Male/Female
ZML-9540	.047"	Male/Male
ZML-9640	.047"	Male/Female

All Zefon Gas Sampling Lines come individually packaged with 50 units per case. Call for custom configurations.

## Attachment Options

Attachments are available pre-installed on any line. To order, select the option you want and adjust the product number to reflect the option. To do this, simply replace the last digit of the product number with the option number. Example: To add option #3 to a ZML-9520, change the number to ZML-9523.

Option #1 - Add a 5.0µm Filter to the sampling line

Option #3 - Add a 12" Nafion® dryer segment to the sampling line

Option #6 - Add both a 5.0µm filter and a 12" Nafion® dryer segment to the line.

Nafion is a registered trademark of Dupont



Distributed By:



www.I-MA.com  
info@I-MA.com  
p: 888.548.4462  
f: 888-548-1462