

StatCheck™



Disposable Blood Pressure Cuffs

For over ten years Zefon International has been designing and manufacturing superior blood pressure cuffs for leading blood pressure measurement instrument manufacturers. With a complete range of sizes, materials and connectors, Zefon brand cuffs represent the best value to bring consistent performance, patient comfort and economy to your medical setting.

Features

- Disposable; Eliminates costly disinfection procedures and contamination concerns.
- Cost Effective; Economical enough for single patient use.
- Increase infection control with single patient use cuffs.
- White vinyl material is easy to wipe off and keep clean.
- Meets AAMI and AHA Cuff Sizing Recommendations.
- Color coded printing on each cuff makes it easy to identify cuff size.
- Latex Free.
- Designed as a direct replacement for OEM cuffs provided by monitor manufacturers.
- Full range of connectors and adapters allow standardization of cuffs across the whole hospital.
- Easy connection to electronic NIBP monitors or manual blood pressure measurement systems.

Latex Free **CE**




5350 SW 1st Lane, Ocala, FL 34474 USA
Tel: +1 3528548080; Fax: +1 3528547480
www.zefon.com

MM00004 Rev 4 2008-06

EC REP Blue Box Medical Limited
Unit 29, New Forest Enterprise Centre,
Rushington Business Park; Chapel Lane
Totton, Southampton, SO40 9LA United Kingdom
Phone: +023 8066 9000; Fax: +023 8066 9002

Ordering Information

Single Tube Cuffs




		Monitor				
		CAS 740/750 Criticare 508 Datascopes Accutor Plus & Passport 2 Philips/HP Siemens Spacelabs	Colin Datascopes Passport XG & Expert Nihon Kodan Lifescope 9 & 12 Welch Allyn Vital Signs/Atlas (Before July 2000)	Welch Allyn Vital Signs/Atlas (after July 2000)	Welch Allyn Protocol Propaq	
Cuff Size	Color Code	Limb Circumference	1-Tube Bayonet Connector	1-Tube Female Luer Connector	1-Tube Screw Type Connector	1-Tube Female Twist Lock Connector
Child	Green	11-22cm	ZSC4401HP	ZSC4401F	ZSC4401D	ZSC4401MQF
Small Adult	Light Blue	16-28cm	ZSC4301HP	ZSC4301F	ZSC4301D	ZSC4301MQF
Adult	Navy Blue	26-33cm	ZSC4201HP	ZSC4201F	ZSC4201D	ZSC4201MQF
Adult Long*	Navy Blue	26-33cm	ZSC4221HP	ZSC4221F	ZSC4221D	ZSC4221MQF
Large Adult	Red	33-41cm	ZSC4101HP	ZSC4101F	ZSC4101D	ZSC4101MQF
Large Adult Long*	Red	33-41cm	ZSC4121HP	ZSC4121F	ZSC4121D	ZSC4121MQF
Thigh	Brown	39-55cm	ZSC4001HP	ZSC4001F	ZSC4001D	ZSC4001MQF

* Available late summer 2008

All cuffs are packaged with 10 cuffs per box

Double Tube Cuffs



		Monitor				
		Critikon Dinamap Criticare 507 Invivo Research Nihon Kodan (Lifescope L) Spacelabs	Datex-Ohmeda	GE Medical Eagle 4000 Tram 100/200/300	Manual	
Cuff Size	Color Code	Limb Circumference	2-Tube Screw Type Connector	2-Tube Male/Female Twist Lock Connector	2-Tube Female Twist Lock Connector	2-Tube Manual Inflation with Male Luer
Child	Green	11-22cm	ZSC4402D	ZSC4402MQMF	ZSC4402MQF	ZSC4412L
Small Adult	Light Blue	16-28cm	ZSC4302D	ZSC4302MQMF	ZSC4302MQF	ZSC4312L
Adult	Navy Blue	26-33cm	ZSC4202D	ZSC4202MQMF	ZSC4202MQF	ZSC4212L
Adult Long*	Navy Blue	26-33cm	ZSC4222D	ZSC4222MQMF	ZSC4222MQF	ZSC4232L
Large Adult	Red	33-41cm	ZSC4102D	ZSC4102MQMF	ZSC4102MQF	ZSC4112L
Large Adult Long*	Red	33-41cm	ZSC4122D	ZSC4122MQMF	ZSC4122MQF	ZSC4132L
Thigh	Brown	39-55cm	ZSC4002D	ZSC4002MQMF	ZSC4002MQF	ZSC4012L

* Available late summer 2008

All cuffs are packaged with 10 cuffs per box

Distributed By:



5350 SW 1st Lane,
Ocala, FL 34474 USA
Tel: +1 3528548080
www.zefon.com



www.I-MA.com
info@I-MA.com
p: 888.548.4462
f: 888-548-1462