

VitalCheck™



Reusable Blood Pressure Cuffs

For over ten years Zefon International has been designing and manufacturing superior blood pressure cuffs for leading blood pressure measurement instrument manufacturers. With a complete range of sizes, materials and connectors, Zefon brand cuffs represent the best value to bring consistent performance, patient comfort and economy to your medical setting.

Features

- Reusable; Durable construction for long life and multiple patient use.
- Cost Effective; Bladder and sleeve design allows sleeves to be replaced as they become worn or contaminated.
- Made from Urethane material that is soft and comfortable for the patient.
- Meets AAMI and AHA Cuff Sizing Recommendations.
- Color coded cuff sleeve makes it easy to identify cuff size.
- Latex Free.
- Designed as a direct replacement for OEM cuffs provided by monitor manufacturers.
- Full range of connectors and adapters allow standardization of cuffs across the whole hospital.
- Easy connection to electronic NIBP monitors or manual blood pressure measurement systems.

Latex Free **CE**




5350 SW 1st Lane, Ocala, FL 34474 USA
Tel: +1 3528548080; Fax: +1 3528547480
www.zefon.com

MM00005 Rev 4 2008-06

EC REP Blue Box Medical Limited
Unit 29, New Forest Enterprise Centre,
Rushington Business Park; Chapel Lane
Totton, Southampton, SO40 9LA United Kingdom
Phone: +023 8066 9000; Fax: +023 8066 9002

Ordering Information


Single Tube Cuffs



Cuff Size	Color Code	Limb Circumference	1-Tube Bayonet Connector	1-Tube Female Luer Connector	1-Tube Screw Type Connector	1-Tube Female Twist Lock Connector
Infant	Brown	9-15cm	ZVC4901HP	ZVC4901F	ZVC4901D	ZVC4901MQF
Child	Green	14-21.5cm	ZVCC4401HP	ZVCC4401F	ZVCC4401D	ZVCC4401MQF
Small Adult	Light Blue	20.5-28.5cm	ZVCC4301HP	ZVCC4301F	ZVCC4301D	ZVCC4301MQF
Adult	Navy Blue	27-35cm	ZVCC4201HP	ZVCC4201F	ZVCC4201D	ZVCC4201MQF
Large Adult	Plum	34-44cm	ZVCC4101HP	ZVCC4101F	ZVCC4101D	ZVCC4101MQF
Thigh	Brown	42-52cm	ZVC4001HP	ZVC4001F	ZVC4001D	ZVC4001MQF

All cuffs are packaged with 5 cuffs per box

Double Tube Cuffs



Cuff Size	Color Code	Limb Circumference	2-Tube Screw Type Connector	2-Tube Male/Female Twist Lock Connector	2-Tube Female Twist Lock Connector	2-Tube Manual Inflation with Male Luer
Infant	Brown	9-15cm	ZVC4902D	ZVC4902MQMF	ZVC4902MQF	ZVC4912L
Child	Green	14-21.5cm	ZVCC4402D	ZVCC4401MQMF	ZVCC4402MQF	ZVCC4412L
Small Adult	Light Blue	20.5-28.5cm	ZVCC4302D	ZVCC4302MQMF	ZVCC4302MQF	ZVCC4312L
Adult	Navy Blue	27-35cm	ZVCC4202D	ZVCC4202MQMF	ZVCC4202MQF	ZVCC4212L
Large Adult	Plum	34-44cm	ZVCC4102D	ZVCC4102MQMF	ZVCC4102MQF	ZVCC4112L
Thigh	Brown	42-52cm	ZVC4002D	ZVC4002MQMF	ZVC4002MQF	ZVC4012L

All cuffs are packaged with 5 cuffs per box

Distributed By:



3530 SW 1st Lane,
Ocala, FL 34474 USA
Tel: +1 3528548080
www.zefon.com



www.I-MA.com
info@I-MA.com
p: 888.548.4462
f: 888-548-1462